


New Colored Flow Moderators



DURECT Corporation introduces a line of colored flow moderators for use with ALZET® Osmotic Pumps during *in vivo* imaging procedures and for pump identification purposes.



In Vivo Imaging

- Eliminate background signal during *in vivo* bio-luminescence imaging using IVIS imaging systems.
- Improve the quality of your *in vivo* imaging data.
- Studies show that colored ALZET Flow Moderators eliminate the high background interference observed from white Flow Moderators during *in vivo* imaging of animals.

Pump Identification

- Easily differentiate between pumps filled with different solutions: control and experimental pumps
- Eliminate the risk of mistakenly implanting a control pump into an animal in the experimental group
- Three colors to choose from: blue, teal, and white[†]

[†] Included with each package of pumps.

Description	Color	Size	Item No.
Color Flow Moderator	Blue	100X	0002609
Color Flow Moderator	Teal	100X	0002607
Color Flow Moderator	Blue	200X	0002489
Color Flow Moderator	Teal	200X	0002488
Color Flow Moderator	Blue	2MLX	0002504
Color Flow Moderator	Teal	2MLX	0002503

[‡] Sold in quantities of 10 per bag

Customer Service
877-922-5938
alzetcs@direct.com

www.alzet.com

Technical Support
800-692-2990
alzet@direct.com